

**Westfalen Medical BV**

Rigastraat 14

7418 EW Deventer










Tel. 0570 - 858 450

Fax 0570 - 858 451

www.westfalenmedical.com

info@westfalenmedical.nl

## Instruction card for oxygen cylinder and regulator

<p>1. Always read the instructions.</p> 	<p>2. Do not smoke. No open fire.</p> 	<p>3. Ventilate the room well.</p> 
<p>4. Ensure that the cylinder is stood upright on a flat surface.</p> 	<p>5. Do not use oil/grease.</p> 	<p>6. Maximum length connection hose.</p> 
<p>7. We recommend replacing the nozzle once a week.</p> 	<p>8. Do not use any electrical equipment.</p> 	<p>9. Malfunctions and/or questions:</p>  <p>Tel. 0570 - 858 450</p>

### Regulator placement

1. Screw the regulator onto the cylinder (hand-tight). Be sure that your hands are clean and free of grease.
2. Open the cylinder by gently turning the main valve of the cylinder. The rotary knob is located on top of the cylinder.
3. Configure the flow meter (number of litres per minute) as prescribed by your physician.
4. In the absence of flow: check the flow by placing the nozzle and/or hose on the palm of your hand or in a glass of water. In the absence of bubbles, inspect the hose and flow meter. If you still do not experience any flow, contact Westfalen Medical.

### Regulator removal

1. Close the cylinder by gently turning the main valve of the cylinder. The rotary knob is located on top of the cylinder.
2. Set the flow meter of the regulator to setting 2 or its highest setting.
3. You will now hear/let residual oxygen flow out of the regulator.
4. Check the contents gauge to see whether all of the oxygen has left the regulator. The gauge should read "0".
5. Set the flow meter to "0".
6. Unscrew the regulator from the cylinder. Use 2 fingers without applying any force. Attention! Do not use any tools for this operation.